


OPTIMUS Clamp System



INTRODUCTION

Overview

Optimus Clamp System, an innovative cranial clamping solution, allows surgeons to do easy, safe and fast fixation of the bone flap to the cranium.

Various clinical challenges surgeons are facing these days in cranial fixation surgeries require solutions that can be only solved with innovative products.

Osteonic proudly meets these high demands, thanks to its cutting-edge development and manufacturing power coming from its key personnel with in-depth knowledge and experiences of more than two decades in products.

- ▶ Clamp type design saves time in operating room by eliminating pre-drilling and screwing process.
- ▶ Various size of clamps(12 mm, 16mm, 20mm)allows surgeons to take various options.
- ▶ Being made of Titanium Alloy, implants free from X-ray and MR.
- ▶ Low profile clamps provide a minimal palpability.
- ▶ Sterile package provides safety and convenience.

Ordering information



12mm OPT-12
OPT-12T(Tilttable)



16mm OPT-16
OPT-16T(Tilttable)



20mm OPT-20
OPT-20T(Tilttable)



Ordering information

Applying Forcep

OPT-001



Holding Forcep

OPT-002



Rod Cutter

OPT-003



Removal Forcep

OPT-004



Wrench

OPT-006



OPT-005





OSTEONIC Co., Ltd.

1206Ho, 38, Digital-ro 29-gil, Guro-gu, Seoul, Korea P.O. #152-779
TEL: +822 6902 8400 FAX: +822 6902 8401 www.osteonic.com

CS - E - 1411